

**PAJUNK®**

## ***Rectus Sheath Kits***


*Catheter System For a  
Continuous Technique  
and Single Shot*

*Truncal Blocks*

*Rectus sheath blocks have been described for many years.*

*Abdominal field blocks since 1899 by Schleich (1) with the development of ultrasound techniques, rectus sheath catheters can now be placed with ultrasound guidance (2,5)*

*Rectus sheath catheters have shown as part of an enhanced recovery/multimodal analgesia technique giving low pain scores and low opiate requirements (3,4)*



*PAJUNK® can now offer varying formats of continuous kits and single shot needles including TuohySono, the ultrasound visible needle.*

### *TuohySono Cannula*

*The TuohySono Cannula from PAJUNK® combines proven Tuohy tip with the unique reflection of the patented “Cornerstone” reflectors. The cannulas are now available in individual packaging.*



## Rectus Sheath Kit

Product	Product Code	Units
<b>Rectus Sheath Kit</b>	0721153-49	10
1 x Touhy needle 16g x 80mm plastic stylet 2 x catheters 19g x 90cm with closed tip 2 x clamping adaptors (including 2 x White Luer Caps) 1 x Y connector (including 3 x White Luer Hub Caps connected to each end of the tubing) 1 x filter 2 x (Rectus Sheath) labels –Yellow- 2 x FixoCath 1 x lancet 1 x Luer Cap White		




## Rectus Sheath Pack

Product	Product Code	Units
<b>Rectus Sheath Pack</b>	071151-927	20
1 x Blunt Drawing Up needle 2 x 20ml Luer slip syringes 2 x 10ml Luer slip Syringes 10 x Green Swabs 1 x fenestrated drape with clear window Outer wrap 1 x Compartment tray		
1 x Touhy needle 16g x 80mm plastic stylet 2 x catheters 19g x 90cm with closed tip 2 x clamping adaptors (including 2 x White Luer Caps) 1 x connector (including 3 x White Luer Hub Caps connected to each end of the tubing) 1 x filter 2 x (Rectus Sheath) labels –Yellow- 2 x FixoCath 1 x lancet 1 x Luer Cap White		

*The TuohySono Cannula has removable wings and is available in three sizes*

Product	Size	Code	Units
TuohySono cannula	17g x 90mm	1085-4M090	25
TuohySono cannula	18g x 90mm	1085-4K090	25
TuohySono cannula	18g x 50mm	1085-4K050	25

1. Schleich CL. *Schmerzlose Operationen*. Berlin: Springer; 1899: 240-1
2. Sandeman DJ, Dilley AV. *Ultrasound-guided rectus sheath block and catheter placement*. *ANZ Journal of Surgery* 2008; 78: 621-3.
3. Cornish P, Deacon A. *Rectus sheath catheters for continuous analgesia after upper abdominal surgery*. *ANZ Journal of Surgery* 2007; 77: 84.
4. Webster, K, Hubble, S. *Rectus sheath analgesia in intensive care patients: technique description and case series*. *Anaesthesia and Intensive Care* 2009; 37: 855.
5. Webster, K. *Ultrasound guided rectus sheath block – analgesia for abdominal surgery*. *Update in Anaesthesia* 2010: 26; 12



PAJUNK GmbH  
 Medizintechnologie  
 Karl-Hall-Strasse 1  
 D-78187 Geisingen/Germany  
 Phone +49 (0) 77 04/92 91-0  
 Fax +49 (0) 77 04/92 91-6 00  
[www.pajunk.com](http://www.pajunk.com)